

Product Datasheet

DYKDDDDK TAG antibody (HRP) (orb436435)

Description	Rat monoclonal antibody against DYKDDDDK TAG conjugated to HRP
Species/Host	Rat
Conjugation	HRP
Tested Applications	ELISA, IHC-P, WB
Target	DYKDDDDK TAG
Preservatives	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Form/Appearance	Purified IgG conjugated to Horseradish Peroxidase (HRP) - liquid
Concentration	1.0mg/ml
Storage	+4°C DO NOT FREEZE
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Expiration Date	12 months from date of receipt.