

Product Datasheet

MOUSE ANTI RABBIT IgM (B CELL MARKER):FITC (orb436468)

Description	Mouse monoclonal antibody against IgM conjugated to FITC
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	FITC
Tested Applications	FC
Target	IgM
Preservatives	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Form/Appearance	Purified IgG conjugated to Fluorescein Isothiocyanate Isomer 1 (FITC) - liquid
Concentration	0.1 mg/ml
Storage	+4°C, -20°C if preferred
Note	For research use only
Application notes	(B-CELL MARKER)
Isotype	IgG1
Clonality	Monoclonal
Clone Number	NRBM
Uniprot ID	P04221, P03988
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)