

Product Datasheet

beta 2 microglobulin antibody (FITC) (orb43656)

Description

Mouse Monoclonal to beta 2 microglobulin.


Reactivity

Human

Conjugation

FITC

Tested Applications

FC

Immunogen

Purified human beta2-microglobulin

Target

beta2-Microglobulin

Preservatives

Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Concentration

1 mg/ml

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: Recommended dilution: 1 µg/ml; positive control: peripheral blood lymphocytes, negative control: DAUDI human Burkitt lymphoma cell line, erythrocytes.

Isotype

Mouse IgG2a

Clonality

Monoclonal

Purity

Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID
[P61769](#)
Entrez
[567](#)
Dilution Range

Flow cytometry: Recommended dilution: 1 µg/ml; positive control: peripheral blood lymphocytes, negative control: DAUDI human Burkitt lymphoma cell line, erythrocytes.

Expiration Date

12 months from date of receipt.

Flow cytometry analysis (surface stainin...