

## Product Datasheet

### Canine Ang1-7 ELISA Kit (orb436843)

<b>Description</b>	Canine Angiotensin 1-7 (Ang1-7) competitive inhibition ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	15.625-1000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Sensitivity</b>	5.33pg/mL
<b>Expiration Date</b>	Please enquire.