

## Product Datasheet

### Canine AngI ELISA Kit (orb436848)

<b>Description</b>	Canine Angiotensin I (AngI) competitive inhibition ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	31.25-2000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
<b>Sensitivity</b>	9.28pg/mL
<b>Expiration Date</b>	Please enquire.