

Product Datasheet

Canine BDNF ELISA Kit (orb436869)

Description	Canine Brain Derived Neurotrophic Factor (BDNF) Sandwich ELISA Kit
Reactivity	Canine
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q7YRB4
Sensitivity	0.055ng/mL
Entrez	403461
Expiration Date	Please enquire.