

## Product Datasheet

### Human C1QTNF1 ELISA Kit (orb436900)

<b>Description</b>	Human C1q And Tumor Necrosis Factor Related Protein 1 (C1QTNF1) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1.56-100ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma and other biological fluids.
<b>Sample Types</b>	serum, plasma and other biological fluids.
<b>Uniprot ID</b>	<a href="#">Q9BXJ1</a>
<b>Sensitivity</b>	0.57ng/mL
<b>Entrez</b>	<a href="#">114897</a>
<b>Expiration Date</b>	Please enquire.