

## Product Datasheet

### Human CSF2Rb ELISA Kit (orb436904)

<b>Description</b>	Human Colony Stimulating Factor 2 Receptor Beta (CSF2Rb) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	0.156-10ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	tissue homogenates, cell lysates and other biological fluids.
<b>Sample Types</b>	tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">P32927</a>
<b>Sensitivity</b>	0.061ng/mL
<b>Entrez</b>	<a href="#">1439</a>
<b>Expiration Date</b>	Please enquire.