

Product Datasheet

Porcine CTGF ELISA Kit (orb436908)

| | |
|--------------------------|--|
| Description | Porcine Connective Tissue Growth Factor (CTGF) Sandwich ELISA Kit |
| Reactivity | Porcine |
| Range | 0.625-40ng/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates and other biological fluids. |
| Uniprot ID | F1S3P6 |
| Sensitivity | 0.242ng/mL |
| Entrez | 100152141 |
| Expiration Date | Please enquire. |