

Product Datasheet

Rabbit CTGF ELISA Kit (orb436910)

Description	Rabbit Connective Tissue Growth Factor (CTGF) Sandwich ELISA Kit
Reactivity	Rabbit
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q2WG38
Sensitivity	18.54pg/mL
Entrez	100008638
Expiration Date	Please enquire.