

Product Datasheet

Bovine EG-VEGF ELISA Kit (orb436931)

Description	Bovine Endocrine Gland Derived Vascular Endothelial Growth Factor (EG-VEGF) Sandwich ELISA Kit
Reactivity	Bovine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	F1MN36
Sensitivity	6.3pg/mL
Entrez	617599
Expiration Date	Please enquire.