

Product Datasheet

Mouse FGF15 ELISA Kit (orb436947)

Description	Mouse Fibroblast Growth Factor 15 (FGF15) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	035622
Sensitivity	12.4pg/mL
Entrez	14170
Expiration Date	Please enquire.