

## Product Datasheet

### Human FGF17 ELISA Kit (orb436949)

<b>Description</b>	Human Fibroblast Growth Factor 17 (FGF17) Sandwich ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	78.125-5000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<b>O60258</b>
<b>Sensitivity</b>	35.96pg/mL
<b>Entrez</b>	<b>8822</b>
<b>Expiration Date</b>	Please enquire.