

Product Datasheet

Porcine FGF2 ELISA Kit (orb436974)

Description	Porcine Fibroblast Growth Factor 2, Basic (FGF2) Sandwich ELISA Kit
Reactivity	Porcine
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	F1RQN3
Sensitivity	41.78pg/mL
Entrez	397643
Expiration Date	Please enquire.