

## Product Datasheet

### Porcine FGF7 ELISA Kit (orb436992)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Porcine Fibroblast Growth Factor 7 (FGF7)<br>Sandwich ELISA Kit                          |
| <b>Reactivity</b>        | Porcine  |
| <b>Range</b>             | 15.625-1000pg/mL   |
| <b>Form/Appearance</b>   | Traditional  |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | serum, plasma, tissue homogenates and other biological fluids.                           |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates and other biological fluids.                           |
| <b>Uniprot ID</b>        | <a href="#">Q9N198</a>   |
| <b>Sensitivity</b>       | 5.9pg/mL   |
| <b>Entrez</b>            | <a href="#">397281</a>   |
| <b>Expiration Date</b>   | Please enquire.  |