

Product Datasheet

Human HBEGF ELISA Kit (orb437026)

| | |
|--------------------------|---|
| Description | Human Heparin Binding Epidermal Growth Factor Like Growth Factor (HBEGF) Sandwich ELISA Kit |
| Reactivity | Human |
| Range | 15.625-1000pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids. |
| Sample Types | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| Uniprot ID | Q99075 |
| Sensitivity | 6.3pg/mL |
| Entrez | 1839 |
| Expiration Date | Please enquire. |