

## Product Datasheet

### Porcine IGFBP3 ELISA Kit (orb437114)

<b>Description</b>	Porcine Insulin Like Growth Factor Binding Protein 3 (IGFBP3) Sandwich ELISA Kit
<b>Reactivity</b>	Porcine
<b>Range</b>	3.125-200ng/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">P16611</a>
<b>Sensitivity</b>	1.485ng/mL
<b>Entrez</b>	<a href="#">448812</a>
<b>Expiration Date</b>	Please enquire.