

## Product Datasheet

### Bovine IL17 ELISA Kit (orb437184)

<b>Description</b>	Bovine Interleukin 17 (IL17) Sandwich ELISA Kit
<b>Reactivity</b>	Bovine
<b>Range</b>	31.25-2000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<b>B3GDZ3</b>
<b>Sensitivity</b>	11.2pg/mL
<b>Entrez</b>	<b>282863</b>
<b>Expiration Date</b>	Please enquire.