

Product Datasheet

Human IL18RAP ELISA Kit (orb437215)

Description	Human Interleukin 18 Receptor Accessory Protein (IL18RAP) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	095256
Sensitivity	0.051ng/mL
Entrez	8807
Expiration Date	Please enquire.