

Product Datasheet

Bovine IL1RA ELISA Kit (orb437249)

Description	Bovine Interleukin 1 Receptor Antagonist (IL1RA) Sandwich ELISA Kit
Reactivity	Bovine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	O77482
Sensitivity	13.5pg/mL
Entrez	281860
Expiration Date	Please enquire.