

Product Datasheet

NCOA6 antibody (Biotin) (orb43730)

Description	NCOA6 antibody (Biotin)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Human Nuclear receptor coactivator 6 protein (285-504AA)
Target	NCOA6
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q14686
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.