

Product Datasheet

Mouse IL5Ra ELISA Kit (orb437348)

Description	Mouse Interleukin 5 Receptor Alpha (IL5Ra) Sandwich ELISA Kit
Reactivity	Mouse
Range	62.5-4000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P21183
Sensitivity	28.5pg/mL
Entrez	16192
Expiration Date	Please enquire.