

Product Datasheet

Guinea pig IL8 ELISA Kit (orb437374)

Description	Guinea pig Interleukin 8 (IL8) Sandwich ELISA Kit
Reactivity	Guinea pig
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P49113
Sensitivity	6.8pg/mL
Entrez	100379599
Expiration Date	Please enquire.