

## Product Datasheet

### Porcine IL8 ELISA Kit (orb437376)

|                          |                                                                                                         |
|--------------------------|---------------------------------------------------------------------------------------------------------|
| <b>Description</b>       | Porcine Interleukin 8 (IL8) Sandwich ELISA Kit                                                          |
| <b>Reactivity</b>        | Porcine                                                                                                 |
| <b>Range</b>             | 31.2-2000pg/mL                                                                                          |
| <b>Form/Appearance</b>   | Traditional                                                                                             |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt.                |
| <b>Note</b>              | For research use only                                                                                   |
| <b>Application notes</b> | serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.   |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids. |
| <b>Uniprot ID</b>        | <b>P26894</b>                                                                                           |
| <b>Sensitivity</b>       | 12.3pg/mL                                                                                               |
| <b>Entrez</b>            | <b>396880</b>                                                                                           |
| <b>Expiration Date</b>   | Please enquire.                                                                                         |

---

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)