

## Product Datasheet

### Mouse IL8Rb ELISA Kit (orb437384)

<b>Description</b>	Mouse Interleukin 8 Receptor Beta (IL8Rb) Sandwich ELISA Kit
<b>Reactivity</b>	Mouse
<b>Range</b>	31.25-2000pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates and other biological fluids.
<b>Uniprot ID</b>	<a href="#">P35343</a>
<b>Sensitivity</b>	12.4pg/mL
<b>Entrez</b>	<a href="#">12765</a>
<b>Expiration Date</b>	Please enquire.