

## Product Datasheet

### Mouse ITaC ELISA Kit (orb437409)

|                          |  |
|--------------------------|--|
| <b>Description</b>       | Mouse Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) Sandwich ELISA Kit        |
| <b>Reactivity</b>        | Mouse  |
| <b>Range</b>             | 31.25-2000pg/mL  |
| <b>Form/Appearance</b>   | Traditional  |
| <b>Storage</b>           | All reagents should be kept according to the instructions on the packaging upon receipt. |
| <b>Note</b>              | For research use only  |
| <b>Application notes</b> | serum, plasma and other biological fluids.   |
| <b>Sample Types</b>      | serum, plasma and other biological fluids.   |
| <b>Uniprot ID</b>        | <a href="#">Q9JHH5</a>   |
| <b>Sensitivity</b>       | 12.4pg/mL  |
| <b>Entrez</b>            | <a href="#">56066</a>  |
| <b>Expiration Date</b>   | Please enquire.  |