

Product Datasheet

Rat ITaC ELISA Kit (orb437410)

Description	Rat Interferon Inducible T-Cell Alpha Chemoattractant (ITaC) Sandwich ELISA Kit
Reactivity	Rat
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q7TNL0
Sensitivity	16.13pg/mL
Entrez	305236
Expiration Date	Please enquire.