

Product Datasheet

Human PDGFAB ELISA Kit (orb437445)

Description	Human Platelet Derived Growth Factor AB (PDGFAB) Sandwich ELISA Kit
Reactivity	Human
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	27pg/mL
Expiration Date	Please enquire.