

Product Datasheet

Porcine PDGFB ELISA Kit (orb437452)

Description	Porcine Platelet Derived Growth Factor Subunit B (PDGFB) Sandwich ELISA Kit
Reactivity	Porcine
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P20034
Sensitivity	13.45pg/mL
Entrez	100126843
Expiration Date	Please enquire.