

Product Datasheet

Human PDGFRL ELISA Kit (orb437461)

Description	Human Platelet Derived Growth Factor Receptor Like Protein (PDGFRL) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q15198
Sensitivity	0.041ng/mL
Entrez	5157
Expiration Date	Please enquire.