

Product Datasheet

Porcine TGFb1 ELISA Kit (orb437488)

Description	Porcine Transforming Growth Factor Beta 1 (TGFb1) Sandwich ELISA Kit
Reactivity	Porcine
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, platelet-poor plasma, tissue homogenates, cell culture supernates and other biological fluids.
Sample Types	serum, platelet-poor plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P07200
Sensitivity	10.5pg/mL
Entrez	397078
Expiration Date	Please enquire.