

Product Datasheet

Mouse TGFb2 ELISA Kit (orb437495)

Description	Mouse Transforming Growth Factor Beta 2 (TGFb2) Sandwich ELISA Kit
Reactivity	Mouse
Range	31.25-2000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell culture supernatants and other biological fluids.
Uniprot ID	P27090
Sensitivity	13.2pg/mL
Entrez	21808
Expiration Date	Please enquire.