

## Product Datasheet

### Canine TNFa ELISA Kit (orb437515)

<b>Description</b>	Canine Tumor Necrosis Factor Alpha (TNFa) Sandwich ELISA Kit
<b>Reactivity</b>	Canine
<b>Range</b>	3.125-200pg/mL
<b>Form/Appearance</b>	Traditional
<b>Storage</b>	All reagents should be kept according to the instructions on the packaging upon receipt.
<b>Note</b>	For research use only
<b>Application notes</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
<b>Sample Types</b>	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
<b>Uniprot ID</b>	<a href="#">P51742</a>
<b>Sensitivity</b>	1.19pg/mL
<b>Entrez</b>	<a href="#">403922</a>
<b>Expiration Date</b>	Please enquire.