

Product Datasheet

Canine TNFb ELISA Kit (orb437529)

Description	Canine Tumor Necrosis Factor Beta (TNFb) Sandwich ELISA Kit
Reactivity	Canine
Range	15.625-1000ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q5WR07
Sensitivity	3.6ng/mL
Entrez	607183
Expiration Date	Please enquire.