

Product Datasheet

Porcine TNFb ELISA Kit (orb437532)

Description	Porcine Tumor Necrosis Factor Beta (TNFb) Sandwich ELISA Kit
Reactivity	Porcine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P26445
Sensitivity	5.1pg/mL
Entrez	407742
Expiration Date	Please enquire.