

Product Datasheet

General 17-OHP ELISA Kit (orb437603)

Description	General 17-Hydroxyprogesterone (17-OHP) competitive inhibition ELISA Kit
Reactivity	All
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Sensitivity	0.092ng/mL
Expiration Date	Please enquire.