

Product Datasheet

General 8-OHdG ELISA Kit (orb437610)

| | |
|--------------------------|--|
| Description | General 8-Hydroxydeoxyguanosine (8-OHdG) competitive inhibition ELISA Kit |
| Reactivity | All |
| Range | 62.5-4000pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma and other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Sensitivity | 18.5pg/mL |
| Expiration Date | Please enquire. |