

## Product Datasheet

### EBV antigen EBNA-1 antibody (orb43764)

<b>Description</b>	Mouse Monoclonal to EBV antigen EBNA-1.
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, IHC-P
<b>Target</b>	EBNA-1 (EBV)
<b>Preservatives</b>	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
<b>Concentration</b>	1 mg/ml
<b>Storage</b>	Store at 2-8°C. Do not freeze.
<b>Note</b>	For research use only
<b>Isotype</b>	Mouse IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	Purified by protein-A affinity chromatography.
<b>Uniprot ID</b>	<a href="#">P03211</a>
<b>Entrez</b>	<a href="#">3783709</a>
<b>Expiration Date</b>	12 months from date of receipt.