

Product Datasheet

Rat ACTH ELISA Kit (orb437693)

| | |
|--------------------------|--|
| Description | Rat Adrenocorticotrophic Hormone (ACTH) competitive inhibition ELISA Kit |
| Reactivity | Rat |
| Range | 15.625-1000pg/mL |
| Form/Appearance | Traditional |
| Storage | All reagents should be kept according to the instructions on the packaging upon receipt. |
| Note | For research use only |
| Application notes | serum, plasma and other biological fluids. |
| Sample Types | serum, plasma and other biological fluids. |
| Sensitivity | 4.76pg/mL |
| Expiration Date | Please enquire. |