

Product Datasheet

Rat ACTHR ELISA Kit (orb437694)

Description	Rat Adrenocorticotrophic Hormone Receptor (ACTHR) Sandwich ELISA Kit
Reactivity	Rat
Range	0.625-40ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	tissue homogenates and other biological fluids.
Sample Types	tissue homogenates and other biological fluids.
Uniprot ID	G3V848
Sensitivity	0.234ng/mL
Entrez	282839
Expiration Date	Please enquire.