

Product Datasheet

Rat ADAMTS1 ELISA Kit (orb437723)

Description	Rat A Disintegrin And Metalloproteinase With Thrombospondin 1 (ADAMTS1) Sandwich ELISA Kit
Reactivity	Rat
Range	78.125-5000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q68EJ2
Sensitivity	30.5pg/mL
Entrez	79252
Expiration Date	Please enquire.