

Product Datasheet

Rat ADAMTS4 ELISA Kit (orb437725)

Description	Rat A Disintegrin And Metalloproteinase With Thrombospondin 4 (ADAMTS4) Sandwich ELISA Kit
Reactivity	Rat
Range	3.125-200ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	Q9ESP7
Sensitivity	1.15ng/mL
Entrez	66015
Expiration Date	Please enquire.