

Product Datasheet

Canine ADM ELISA Kit (orb437751)

Description	Canine Adrenomedullin (ADM) competitive inhibition ELISA Kit
Reactivity	Canine
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	O77559
Sensitivity	4.74pg/mL
Entrez	403817
Expiration Date	Please enquire.