

Product Datasheet

Mouse ADM ELISA Kit (orb437753)

Description	Mouse Adrenomedullin (ADM) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P97297
Sensitivity	5.31pg/mL
Entrez	11535
Expiration Date	Please enquire.