

Product Datasheet

Human AGR2 ELISA Kit (orb437799)

Description	Human Anterior Gradient Protein 2 (AGR2) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	O95994
Sensitivity	0.061ng/mL
Entrez	10551
Expiration Date	Please enquire.