

Product Datasheet

Rat ALB ELISA Kit (orb437838)

Description	Rat Albumin (ALB) competitive inhibition ELISA Kit
Reactivity	Rat
Range	3.12-200µg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	P02770
Sensitivity	0.93µg/mL
Entrez	24186
Expiration Date	Please enquire.