

Product Datasheet

Human AMH ELISA Kit (orb437888)

Description	Human Anti-Mullerian Hormone (AMH) competitive inhibition ELISA Kit
Reactivity	Human
Range	46.875-3000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	P03971
Sensitivity	23.2pg/mL
Entrez	268
Expiration Date	Please enquire.