

Product Datasheet

Mouse AMH ELISA Kit (orb437889)

Description	Mouse Anti-Mullerian Hormone (AMH) competitive inhibition ELISA Kit
Reactivity	Mouse
Range	0.312-20ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates, cell lysates, cell culture supernates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernatants and other biological fluids.
Uniprot ID	Q5EC55
Sensitivity	0.085ng/mL
Entrez	11705
Expiration Date	Please enquire.