

Product Datasheet

Human ANXA3 ELISA Kit (orb437951)

Description	Human Annexin A3 (ANXA3) Sandwich ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma, tissue homogenates and other biological fluids.
Sample Types	serum, plasma, tissue homogenates and other biological fluids.
Uniprot ID	P12429
Sensitivity	0.053ng/mL
Entrez	306
Expiration Date	Please enquire.