

Product Datasheet

Human APLN ELISA Kit (orb437981)

Description	Human Apelin (APLN) competitive inhibition ELISA Kit
Reactivity	Human
Range	15.625-1000pg/mL
Form/Appearance	Traditional
Storage	All reagents should be kept according to the instructions on the packaging upon receipt.
Note	For research use only
Application notes	serum, plasma and other biological fluids.
Sample Types	serum, plasma and other biological fluids.
Uniprot ID	Q9ULZ1
Sensitivity	6.13pg/mL
Entrez	8862
Expiration Date	Please enquire.